

Human Varicella-Zoster Virus IgG ELISA Kit

Human Varicella-Zoster Virus IgG ELISA Kit

CATALOG #
AAA204608

REACTIVITY
Human

PRODUCT SYNONYMS

Varicella-Zoster Virus IgG; Human Varicella-Zoster Virus IgG ELISA Kit; Varicella-Zoster Virus IgG elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	12.5ng/ml-200ng/ml
Sensitivity	1 ng/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Varicella-Zoster Virus IgG elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration ofVARICELLA-ZOSTER VIRUS IGGin the sample, This VARICELLA-ZOSTER VIRUS IGGELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versusVARICELLA-ZOSTER VIRUS IGGconcentration. The concentration ofVARICELLA-ZOSTER VIRUS IGGin the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.