

Human Serine/Threonine-Protein Kinase DCLK1 ELISA Kit | DCLK1 elisa kit

Human Serine/Threonine-Protein Kinase DCLK1 ELISA Kit

CATALOG #
AAA205035

REACTIVITY
Human


PRODUCT SYNONYMS

Serine/Threonine-Protein Kinase DCLK1; Human Serine/Threonine-Protein Kinase DCLK1 ELISA Kit; DCLK1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	422
Sensitivity	1 pg/mL
Intra-assay Precision	Intra-assay CV is less than 15%.
Inter-assay Precision	Inter-assay CV is less than 15%.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for DCLK1 elisa kit

Intended Uses: This SERINE/THREONINE-PROTEIN KINASE DCLK1 ELISA kit is intended Laboratory for Research use only. The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of SERINE/THREONINE-PROTEIN KINASE DCLK1 in the sample, this SERINE/THREONINE-PROTEIN KINASE DCLK1 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus SERINE/THREONINE-PROTEIN KINASE DCLK1 concentration. The concentration of SERINE /THREONINE - PROTEIN KINASE DCLK1 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	306518604
NCBI Accession #	NP_001182344.1

NCBI GenBank Nucleotide #	NP_001182344.1
UniProt Accession #	O15075
NCBI Official Full Name	serine/threonine-protein kinase DCLK1 isoform 2
UniProt Protein Name	Serine/threonine-protein kinase DCLK1
UniProt Gene Name	DCLK1
UniProt Synonym Gene Names	DCAMKL1; DCDC3A; KIAA0369
UniProt Entry Name	DCLK1_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/serine-threonine-protein-kinase-dclk1/205035>