

## Human High Density Lipoprotein Serum Amyloid A ELISA Kit

Human High Density Lipoprotein Serum Amyloid A ELISA Kit

CATALOG #  
AAA205659

REACTIVITY  
Human


### PRODUCT SYNONYMS

High Density Lipoprotein Serum Amyloid A; Human High Density Lipoprotein Serum Amyloid A ELISA Kit; High Density Lipoprotein Serum Amyloid A elisa kit

### PRODUCT SPECIFICATIONS

|                 |                                                                   |
|-----------------|-------------------------------------------------------------------|
| Reactivity      | Human                                                             |
| Samples         | Serum, plasma, cell culture supernate and other biological fluids |
| Assay Type      | Sandwich                                                          |
| Detection Range | 1.25 mug/ml-20 mug/ml                                             |
| Sensitivity     | 0.1 mug/ml                                                        |

### IMAGES



Standard Curve (Sample)

### ADDITIONAL INFORMATION

#### Related Product Information for High Density Lipoprotein Serum Amyloid A elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of HIGH DENSITY LIPOPROTEIN SERUM AMYLOID A in the sample, This HIGH DENSITY LIPOPROTEIN SERUM AMYLOID A ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus HIGH DENSITY LIPOPROTEIN SERUM AMYLOID A concentration. The concentration of HIGH DENSITY LIPOPROTEIN SERUM AMYLOID A in the samples is then determined by comparing the O.D. of the samples to the standard curve.

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)

Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/high-density-lipoprotein-serum-amyloid-a/205659>