

Human Soluble E-Cadherin ELISA Kit

Human Soluble E-Cadherin ELISA Kit

CATALOG #
AAA205825

REACTIVITY
Human


PRODUCT SYNONYMS

Soluble E-Cadherin; Human Soluble E-Cadherin ELISA Kit; Soluble E-Cadherin elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sensitivity	1.0 ng/mL
Intra-assay Precision	Intra-assay CV are less than 15%.
Inter-assay Precision	Inter-assay CV are less than 15%.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Soluble E-Cadherin elisa kit

Intended Uses: This sE-Cad ELISA kit is intended for Laboratory for Research use only and is not for use in test or therapeutic procedures. The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of sE-Cad in the sample, this sE-Cad ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus sE-Cad concentration. The concentration of sE-Cad in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.