

Rat Ceramide ELISA Kit

Rat Ceramide ELISA Kit

CATALOG #
AAA209893

REACTIVITY
Rat


PRODUCT SYNONYMS

Ceramide; Rat Ceramide ELISA Kit; Ceramide elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rat
Sensitivity	0.1umol/L
Intra-assay Precision	Intra-assay CV is less than 15%.
Inter-assay Precision	Inter-assay CV is less than 15%.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Ceramide elisa kit

Intended Uses: This CERAMIDE ELISA kit is intended Laboratory for Research use only and is not for use in test or therapeutic procedures. The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of CERAMIDE in the sample, this CERAMIDE ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus CERAMIDE concentration. The concentration of CERAMIDE in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/rat-elisa-kits/ceramide/209893>