

## Rat Tissue Plamnipen Activator ELISA Kit

Rat Tissue Plamnipen Activator ELISA Kit

CATALOG #  
**AAA210002**

REACTIVITY  
**Rat**

### PRODUCT SYNONYMS

Tissue Plamnipen Activator; Rat Tissue Plamnipen Activator ELISA Kit; Tissue Plamnipen Activator elisa kit

### PRODUCT SPECIFICATIONS

|                 |                                                                   |
|-----------------|-------------------------------------------------------------------|
| Reactivity      | Rat                                                               |
| Sequence Length | 166                                                               |
| Samples         | Serum, plasma, cell culture supernate and other biological fluids |
| Assay Type      | Sandwich                                                          |
| Detection Range | 1.25ng/ml-20ng/ml                                                 |
| Sensitivity     | 0.1 ng/ml                                                         |

### IMAGES



Standard Curve (Sample)

### ADDITIONAL INFORMATION

#### Related Product Information for Tissue Plamnipen Activator elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of TISSUE PLAMNIPEN ACTIVATOR in the sample, This TISSUE PLAMNIPEN ACTIVATOR ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus TISSUE PLAMNIPEN ACTIVATOR concentration. The concentration of TISSUE PLAMNIPEN ACTIVATOR in the samples is then determined by comparing the O.D. of the samples to the standard curve.

### NCBI AND UNIPROT INFORMATION

|                         |                                       |
|-------------------------|---------------------------------------|
| NCBI GI #               | 1515452                               |
| NCBI Official Full Name | tissue plasminogen activator, partial |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://www.aaabiotech.com/rat-elisa-kits/tissue-plamnipen-activator/210002>