

Mouse Interferon gamma (IFN-gamma) ELISA Kit | IFN-gamma elisa kit

Mouse Interferon gamma (IFN-gamma) quantitative detection ELISA Kit

CATALOG #
AAA210634

REACTIVITY
Mouse


PRODUCT SYNONYMS

Interferon gamma (IFN-gamma); Mouse Interferon gamma (IFN-gamma) quantitative detection ELISA Kit; Mouse Interferon gamma (IFN-gamma); IFN-gamma elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	This kit identifies natural and recombinant Mouse IFN-gamma, with No crossing with structural analogs.
Assay Type	Quantitative Sandwich
Samples	Serum, plasma, and cell culture supernatant
Detection Range	25 pg/mL - 800 pg/mL
Sensitivity	1.0 pg/ml
Intra-assay Precision	In the same lot CV%<10% 3 groups of known high, medium, and low concentration samples were evaluated 20 times on the same plate.
Inter-assay Precision	Different lot CV%<15: 3 groups of known high, medium, and low concentration samples were evaluated 20 times on different plates.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for IFN-gamma elisa kit	Principle of the Assay: This ELISA kit applies the two steps quantitative sandwich enzyme immunoassay technique. The microtiter plate has been pre-coated with mouse IFN-gamma antibody (solid phase antigen), standards or samples are added to the microtiter plate wells, and the microtiter plate undergoes incubation. Add biotin antibody, incubate, wash thoroughly, and then add avidin enzyme. Wells are thoroughly washed to remove all unbound components, and a sandwich complex of solid antigen-antibody enzyme-labeled secondary antibodies is formed on the solid surface of the microplate. Next, Substrate solution A and substrate solution B are added to each well, and it will produce blue materials under HRP catalysis. The enzyme-substrate reaction is terminated by the addition of sulphuric acid and the color change to yellow. Measured spectrophotometrically at a wavelength of 450 nm, the absorbance (OD value) is positively correlated with the concentration of mouse IFN-gamma in the sample. The concentration of IFN-gamma in the samples could be determined by comparing the O.D. of the samples to the standard curve. Intended Uses: This kit is used for detecting the Mouse IFN-gamma concentration of serum, plasma, and cell culture supernatant.
---	---

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/mouse-elisa-kits/interferon-gamma-ifn-gamma/210634>