

## Pemafibrate Agonist

Pemafibrate

CATALOG #  
**AAA210707**


### PRODUCT SYNONYMS

Pemafibrate; Pemafibrate agonist

### PRODUCT SPECIFICATIONS

|                         |                                                                                    |
|-------------------------|------------------------------------------------------------------------------------|
| Purity/Purification     | >98%                                                                               |
| Solubility              | DMSO                                                                               |
| CAS No.                 | 848259-27-8                                                                        |
| Formula                 | C <sub>28</sub> H <sub>30</sub> N <sub>2</sub> O <sub>6</sub>                      |
| SMILES                  | <chem>CC[C@@H](OC1=CC=CC(CN(C2=NC=CC=CC=C3O2)CCCOC4=CC=C(OC)C=C4)=C1)C(O)=O</chem> |
| Availability            | We cannot ship this item to Switzerland, Italy, Germany, France and Austria.       |
| Preparation and Storage | Store at -20 degree C for ~2 years                                                 |

### IMAGES



Chemical Structure

### ADDITIONAL INFORMATION

|                                                     |                                                                  |
|-----------------------------------------------------|------------------------------------------------------------------|
| Related Product Information for Pemafibrate agonist | Pemafibrate, also known as (R)-K 13675, is a PPAR alpha agonist. |
| Product Categories/Family for Pemafibrate agonist   | <a href="#">PPAR alpha agonist</a>                               |

### NCBI AND UNIPROT INFORMATION

|                  |        |
|------------------|--------|
| Molecular Weight | 490.56 |
|------------------|--------|

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.