

Human SNAP Associated Protein (SNAPAP) RTU ELISA Kit | SNAPAP rtu elisa kit

Human SNAP Associated Protein (SNAPAP) ELISA Kit

CATALOG #
AAA218492

REACTIVITY
Human


PRODUCT SYNONYMS

SNAP Associated Protein (SNAPAP); Human SNAP Associated Protein (SNAPAP) ELISA Kit; SNAPIN; SNAP-25-Binding Protein; Biogenesis of lysosome-related organelles complex 1 subunit 7; Synaptosomal-associated protein 25-binding protein; SNAPAP rtu elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	136
Assay Type	Sandwich
Samples	Tissue Homogenates or Other Biological Fluids.
Detection Range	0.156-10ng/mL
Sensitivity	0.066ng/mL
RTU ELISA Kit's Advantage	The Ready-To-Use (RTU) Kits require fewer steps in setup and take less time to run. Detection Solutions are pre-diluted and total incubation time is shorter than many other ELISA Kits.
Preparation and Storage	For unopened kits: All the reagents should be kept at 4 degree C upon receipt. For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for SNAPAP rtu elisa kit

Test Principle: The microtiter plate provided in this kit has been pre-coated with an antibody specific to SNAPAP. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to SNAPAP. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain SNAPAP, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of SNAPAP in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	6912674
NCBI GeneID	23557
NCBI Accession #	NP_036569.1
NCBI GenBank Nucleotide #	NM_012437.5
UniProt Accession #	O95295; Q5SXU8; D3DV56
Molecular Weight	14,874 Da
NCBI Official Full Name	SNARE-associated protein Snapin
NCBI Official Synonym Full Names	SNAP associated protein
NCBI Official Symbol	SNAPIN
NCBI Official Synonym Symbols	BLOS7; BORCS3; SNAPAP; BLOC1S7
NCBI Protein Information	SNARE-associated protein Snapin
UniProt Protein Name	SNARE-associated protein Snapin
UniProt Gene Name	SNAPIN
UniProt Synonym Gene Names	BLOC1S7; SNAP25BP; SNAPAP; BLOC-1 subunit 7; SNAP-associated protein

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/human-rtu-elisa-kits/snap-associated-protein-snapap/218492>