

Mouse Integrin Associated Protein (IAP) ELISA Kit | Akp3 elisa kit

Mouse Integrin Associated Protein (IAP) ELISA Kit

CATALOG #
AAA219797

REACTIVITY
Mouse


PRODUCT SYNONYMS

Integrin Associated Protein (IAP); Mouse Integrin Associated Protein (IAP) ELISA Kit; CD47; MER6; OA3; Antigen Identified By Monoclonal Antibody 1D8; Antigenic Surface Determinant Protein OA3; Rh-Related Antigen; Akp3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	This assay has high sensitivity and excellent specificity for detection of IAP. No significant cross-reactivity or interference between IAP and analogues was observed.
Assay Type	Quantitative Sandwich
Samples	Mouse tissue homogenates, cell lysates and other biological fluids
Detection Range	0.312-20ng/mL
Sensitivity	< 0.118ng/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level IAP were tested 20 times on one plate, respectively. Intra-Assay: CV<10%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level IAP were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	Shelf Life: 12 months

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Akp3 elisa kit	Principle of the Assay: The microtiter plate provided in this kit has been pre-coated with an antibody specific to IAP. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to IAP. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain IAP, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of IAP in the samples is then determined by comparing the O.D. of the samples to the standard curve. Intended Uses: The kit is a sandwich enzyme immunoassay for the in vitro quantitative measurement of IAP in mouse tissue homogenates, cell lysates and other biological fluids.
Product Categories/Family for Akp3 elisa kit	Traditional ELISA Kit

NCBI AND UNIPROT INFORMATION

NCBI GI #	130745
NCBI GeneID	11648
UniProt Accession #	P24822
Molecular Weight	60,285 Da
NCBI Official Full Name	Intestinal-type alkaline phosphatase
NCBI Official Synonym Full Names	alkaline phosphatase 3, intestine, not Mn requiring
NCBI Official Symbol	Akp3
NCBI Official Synonym Symbols	IAP; dIAP; Akp-3
NCBI Protein Information	intestinal-type alkaline phosphatase
UniProt Protein Name	Intestinal-type alkaline phosphatase
UniProt Gene Name	lap
UniProt Synonym Gene Names	Akp-3; Akp3; IAP; Intestinal alkaline phosphatase

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/mouse-elisa-kits/integrin-associated-protein-iap/219797>