

Antimicrobial Bronidox Accessory Reagent

Antimicrobial Bronidox solution

CATALOG #
AAA224679

PRODUCT SYNONYMS

Antimicrobial Bronidox; Antimicrobial Bronidox solution; 5-Bromo-5-nitro-1,3-dioxane (solution); Bronidox; Antimicrobial Bronidox accessory reagent

PRODUCT SPECIFICATIONS

Form/Format	Transparent non-viscous liquid, 0.2 um filtered, pH 6.0 with 5-Bromo-5-nitro-1,3-dioxane, 5000 ppm, in deionized water.
Concentration	Supplied as a 0.5% aqueous solution (varies by lot)
Additional Remarks	Example of making V ml buffer preserved with e.g. 0.05% BND: Mix 1/10 V of BND with 9/10 V deionized water. Add salts and/or other buffering reagents as per given buffer formula, mix, measure and adjust pH. Filter through 0.2 um filter unit.
Preparation and Storage	Ships on blue ice. Upon receipt, store at room temperature. NOTE: At temperatures below 18°C, product may crystallize out of the saturated solution. DO NOT STORE below 18°C.

ADDITIONAL INFORMATION

Related Product Information for Antimicrobial Bronidox accessory reagent	5-Bromo-5-nitro-1,3-dioxane is an antimicrobial chemical compound and causes inhibition of enzyme activity in bacteria.
Product Categories/Family for Antimicrobial Bronidox accessory reagent	Drugs & Toxicology ; Mounting Media & Preservatives

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/accessory-reagent/antimicrobial-bronidox/224679>