

Buprenorphine Hapten Conjugate

Conjugated Buprenorphine-BSA hapten

CATALOG #
AAA224682

APPLICATIONS
Lateral Flow, ELISA

PRODUCT SYNONYMS

Buprenorphine; Conjugated Buprenorphine-BSA hapten; Buprenorphine-BSA; Buprenex-BSA; Buprenorphine-BSA conjugate; Buprenorphine-BSA hapten; Buprenorphine hapten conjugate

PRODUCT SPECIFICATIONS

Purity/Purification	> 95% pure
Form/Format	Chromatography purified liquid in 10mM PBS, pH 7.4 with 0.1% NaN ₃
Concentration	12.69mg/ml (varies by lot)
Applicable Applications for Buprenorphine hapten conjugate	LF (Lateral Flow), ELISA
Protein Type	Synthetic
Tag/Conjugate	BSA
Assay Information	Matched pair antibody available for Buprenorphine-BSA:
Preparation and Storage	Short Term: 4°C. Long Term: -20°C. Avoid freeze/thaw cycles.

ADDITIONAL INFORMATION

Related Product Information for Buprenorphine hapten conjugate	Buprenorphine is a semi-synthetic opioid that is used to treat opioid addiction in higher dosages (>2 mg), to control moderate acute pain in non-opioid-tolerant individuals in lower dosages (~200 ug), and to control moderate chronic pain in dosages ranging from 20–70 ug/hour. It is available in a variety of formulations.
Product Categories/Family for Buprenorphine hapten conjugate	Drugs & Toxicology ; Diluents ; Synthetic

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/hapten-conjugate/buprenorphine/224682>