

## Histamine Hapten Conjugate

BSA conjugated Histamine Hapten

CATALOG #  
**AAA224688**

### PRODUCT SYNONYMS

Histamine; BSA conjugated Histamine Hapten; Histamine-BSA antigen; Histamine-BSA; Hapten conjugate; Histamine BSA; Histaminum hapten; HIST hapten; Histamine hapten conjugate

### PRODUCT SPECIFICATIONS

|                         |                                                                                                                                                          |
|-------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|
| Purity/Purification     | Purification Method: Dialysis                                                                                                                            |
| Form/Format             | Synthetic histamine conjugated to BSA, purified by dialysis, supplied lyophilized from 0.01 mol/L PBS, pH 7.4.<br>Reconstitute in 0.01 mol/L pbs, pH 7.4 |
| Concentration           | 1.0mg, lyophilized (varies by lot)                                                                                                                       |
| Preparation and Storage | Upon receipt, store at 4°C short term or -20°C long term. Avoid repeat freeze-thaws.                                                                     |

### ADDITIONAL INFORMATION

|                                                             |                                          |
|-------------------------------------------------------------|------------------------------------------|
| Product Categories/Family for<br>Histamine hapten conjugate | Signal Transduction; Diluents; Synthetic |
|-------------------------------------------------------------|------------------------------------------|

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://www.aaabiotech.com/hapten-conjugate/histamine/224688>