

## Galbulin Biochemical

Galbulin

CATALOG #  
**AAA226896**


### PRODUCT SYNONYMS

Galbulin; Galbulin biochemical


### PRODUCT SPECIFICATIONS

|                         |                                                                          |
|-------------------------|--------------------------------------------------------------------------|
| Purity/Purification     | 98%                                                                      |
| Solubility              | <1mg/ml refers to the product slightly soluble or insoluble              |
| Preparation and Storage | Powder: -20 degree C for 3 years<br>In solvent: -80 degree C for 2 years |

### IMAGES



Structure  
CU-76, CAS 2400954-58-5



Structure  
CU-76, CAS 2400954-58-5

### ADDITIONAL INFORMATION

|            |                                                                                                                                                                                                             |
|------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| References | <p>ef="javascript:void(0)" &gt;</p> <p>1. Padilla-Salinas, R., Sun, L., Anderson, R., et al. Discovery of small-molecule cyclic GMP-AMP synthase inhibitors. <i>J. Org. Chem.</i> 85(3):1579-1600(2020)</p> |
|------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

### NCBI AND UNIPROT INFORMATION

|                  |       |
|------------------|-------|
| Molecular Weight | 407.1 |
|------------------|-------|

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.