

## Pomalidomide-amido-C3-COOH Biochemical

Pomalidomide-amido-C3-COOH

CATALOG #  
**AAA228181**


### PRODUCT SYNONYMS

Pomalidomide-amido-C3-COOH; Pomalidomide-amido-C3-COOH biochemical


### PRODUCT SPECIFICATIONS

|                         |                                                                          |
|-------------------------|--------------------------------------------------------------------------|
| Purity/Purification     | 98%                                                                      |
| CAS                     | 2162120-77-4                                                             |
| Solubility              | <1mg/ml refers to the product slightly soluble or insoluble              |
| Preparation and Storage | Powder: -20 degree C for 3 years<br>In solvent: -80 degree C for 2 years |

### IMAGES



Structure  
MPO-IN-1, CAS 2471981-21-0



Structure  
MPO-IN-1, CAS 2471981-21-0

### ADDITIONAL INFORMATION

|                                                                        |                                                                                                                                                                                |
|------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Related Product Information for Pomalidomide-amido-C3-COOH biochemical | Pomalidomide-amido-C3-COOH is a synthesized E3 ligase ligand-linker conjugate that incorporates the Pomalidomide based cereblon ligand and a linker used in PROTAC technology. |
| References                                                             | 1. Patnaik A, et al. Discovery of a novel indole pharmacophore for the irreversible inhibition of myeloperoxidase (MPO). Bioorg Med Chem. 2020;28(12):115548.                  |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.