

## Human Staphylococcus Aureus Enterotoxins (SE) ELISA Kit | SE elisa kit

Human Staphylococcus Aureus Enterotoxins (SE) ELISA Kit

CATALOG #  
**AAA23629**

REACTIVITY  
**Human**


### PRODUCT SYNONYMS

Staphylococcus Aureus Enterotoxins (SE); Human Staphylococcus Aureus Enterotoxins (SE) ELISA Kit; SE elisa kit

### PRODUCT SPECIFICATIONS

|                       |                                  |
|-----------------------|----------------------------------|
| Reactivity            | Human                            |
| Sensitivity           | 0.1 ng/mL                        |
| Intra-assay Precision | Intra-assay CV is less than 15%. |
| Inter-assay Precision | Inter-assay CV is less than 15%. |

### IMAGES



Standard Curve (Sample)

### ADDITIONAL INFORMATION

#### Related Product Information for SE elisa kit

**Intended Uses:** This SE ELISA kit is intended for Laboratory for Research use only and is not for use in diagnostic or therapeutic procedures. In order to determine whether contains SE in the sample, this SE ELISA Kit includes Negative control and Positive control. The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. The color depth was positively correlated with the SE in the sample.

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)

Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/staphylococcus-aureus-enterotoxins-se/23629>