

Human N acetylmuramoyl L alanine amidase (PGLYRP2) ELISA Kit | PGLYRP2 elisa kit

Human N acetylmuramoyl L alanine amidase (PGLYRP2) ELISA Kit

CATALOG #
AAA244967

REACTIVITY
Human


PRODUCT SYNONYMS

N acetylmuramoyl L alanine amidase (PGLYRP2); Human N acetylmuramoyl L alanine amidase (PGLYRP2) ELISA Kit; PGLYRP2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of PGLYRP2. No significant cross-reactivity or interference between PGLYRP2 and analogues was observed. NOTE: Limited by current skills and knowledge, it is impossible for us to complete the cross-reactivity detection between PGLYRP2 and all the analogues, therefore, cross reaction may still exist in some cases.
Samples	Serum, plasma, cell culture supernatants, body fluid and tissue homogenate
Assay Type	Quantitative Competitive
Sensitivity	0.1 ng/mL
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PGLYRP2 elisa kit	Intended Uses: This PGLYRP2 ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Human PGLYRP2. This ELISA kit for research use only, not for therapeutic or test applications! Principle of the Assay: PGLYRP2 ELISA kit applies the competitive enzyme immunoassay technique utilizing a polyclonal anti-PGLYRP2 antibody and an PGLYRP2-HRP conjugate. The assay sample and buffer are incubated together with PGLYRP2-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the PGLYRP2 concentration since PGLYRP2 from samples and PGLYRP2-HRP conjugate compete for the anti-PGLYRP2 antibody binding site. Since the number of sites is limited, as more sites are occupied by PGLYRP2 from the sample, fewer sites are left to bind PGLYRP2-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The PGLYRP2 concentration in each sample is interpolated from this standard curve.
Product Categories/Family for PGLYRP2 elisa kit	Immunology

NCBI AND UNIPROT INFORMATION

NCBI GI #	332982048
NCBI Accession #	YP_004463489.1
NCBI GenBank Nucleotide #	NC_015520.1
Molecular Weight	26,019 Da
NCBI Official Full Name	N-acetylmuramoyl-L-alanine amidase CwID
UniProt Protein Name	N-acetylmuramoyl-L-alanine amidase CwID
UniProt Entry Name	F3ZXZ7_MAHA5

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/n-acetylmuramoyl-l-alanine-amidase-pglyrp2/244967>