

Canine CASP8 (Caspase 8) ELISA Kit | CASP8 elisa kit

Canine CASP8 (Caspase 8) ELISA Kit

CATALOG #
AAA248155

REACTIVITY
Canine


PRODUCT SYNONYMS

CASP8 (Caspase 8); Canine CASP8 (Caspase 8) ELISA Kit; CASP8 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Sequence Length	479
Assay Type	Sandwich ELISA, Double Antibody
Samples	Serum, Plasma, Cell Culture Supernatant, cell or tissue lysate, Other liquid samples
Detection Range	0.156-10ng/ml
Sensitivity	0.094ng/ml
Preparation and Storage	Store entire kit at 2-8C for short-term. For longer-term, please store the microplate & standard at -20C, while the remaining reagents can be stored at 2-8C

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CASP8 elisa kit

Principle of the Assay: This kit was based on sandwich enzyme-linked immune-sorbent assay technology. Anti-CASP8 antibody was pre-coated onto 96-well plates. And the biotin conjugated anti-CASP8 antibody was used as detection antibodies. The standards, test samples and biotin conjugated detection antibody were added to the wells subsequently, and wash with wash buffer. HRP-Streptavidin was added and unbound conjugates were washed away with wash buffer. TMB substrates were used to visualize HRP enzymatic reaction. TMB was catalyzed by HRP to produce a blue color product that changed into yellow after adding acidic stop solution. The density of yellow is proportional to the CASP8 amount of sample captured in plate. Read the O.D. absorbance at 450nm in a microplate reader, and then the concentration of CASP8 can be calculated.

NCBI AND UNIPROT INFORMATION

NCBI GI #

4583149

NCBI GeneID	841
Molecular Weight	61,836 Da
NCBI Official Full Name	caspase 8
NCBI Official Synonym Full Names	caspase 8
NCBI Official Symbol	CASP8
NCBI Official Synonym Symbols	CAP4; MACH; MCH5; FLICE; ALPS2B; Casp-8
NCBI Protein Information	caspase-8
UniProt Protein Name	Caspase-8
UniProt Gene Name	CASP8
UniProt Synonym Gene Names	MCH5; CASP-8; FLICE; MACH

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/canine-elisa-kits/casp8-caspase-8/248155>