

Human D2HGDH (D2-Hydroxyglutarate Dehydrogenase) ELISA Kit | D2HGDH elisa kit

Human D2HGDH (D2-Hydroxyglutarate Dehydrogenase) ELISA Kit

CATALOG #
AAA248174

REACTIVITY
Human


PRODUCT SYNONYMS

D2HGDH (D2-Hydroxyglutarate Dehydrogenase); Human D2HGDH (D2-Hydroxyglutarate Dehydrogenase) ELISA Kit; D2-Hydroxyglutarate Dehydrogenase; D2HGDH elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	410
Assay Type	Sandwich ELISA, Double Antibody
Samples	Serum, Plasma, Cell Culture Supernatant, cell or tissue lysate, Other liquid samples
Detection Range	0.156-10ng/ml
Sensitivity	0.094ng/ml
Preparation and Storage	Store entire kit at 2-8C for short-term. For longer-term, please store the microplate & standard at -20C, while the remaining reagents can be stored at 2-8C

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for D2HGDH elisa kit

Principle of the Assay: This kit was based on sandwich enzyme-linked immune-sorbent assay technology. Anti-D2HGDH antibody was pre-coated onto 96-well plates. And the biotin conjugated anti-D2HGDH antibody was used as detection antibodies. The standards, test samples and biotin conjugated detection antibody were added to the wells subsequently, and wash with wash buffer. HRP-Streptavidin was added and unbound conjugates were washed away with wash buffer. TMB substrates were used to visualize HRP enzymatic reaction. TMB was catalyzed by HRP to produce a blue color product that changed into yellow after adding acidic stop solution. The density of yellow is proportional to the D2HGDH amount of sample captured in plate. Read the O.D. absorbance at 450nm in a microplate reader, and then the concentration of D2HGDH can be calculated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	1036277043
NCBI Official Full Name	D-3-phosphoglycerate dehydrogenase

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/d2hgdh-d2-hydroxyglutarate-dehydrogenase/248174>