

Canine BALP (Bone Alkaline Phosphatase) ELISA Kit | BALP elisa kit

Canine BALP (Bone Alkaline Phosphatase) ELISA Kit

CATALOG #
AAA248344

REACTIVITY
Canine


PRODUCT SYNONYMS

BALP (Bone Alkaline Phosphatase); Canine BALP (Bone Alkaline Phosphatase) ELISA Kit; BALP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Assay Type	Sandwich ELISA, Double Antibody
Samples	Serum, Plasma, Cell Culture Supernatant, cell or tissue lysate, Other liquid samples
Detection Range	0.469-30ng/ml
Sensitivity	0.281ng/ml
Preparation and Storage	Store entire kit at 2-8C for short-term. For longer-term, please store the microplate & standard at -20C, while the remaining reagents can be stored at 2-8C

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for BALP elisa kit

Principle of the Assay: This kit was based on sandwich enzyme-linked immune-sorbent assay technology. Anti-BALP antibody was pre-coated onto 96-well plates. And the biotin conjugated anti-BALP antibody was used as detection antibodies. The standards, test samples and biotin conjugated detection antibody were added to the wells subsequently, and wash with wash buffer. HRP-Streptavidin was added and unbound conjugates were washed away with wash buffer. TMB substrates were used to visualize HRP enzymatic reaction. TMB was catalyzed by HRP to produce a blue color product that changed into yellow after adding acidic stop solution. The density of yellow is proportional to the BALP amount of sample captured in plate. Read the O.D. absorbance at 450nm in a microplate reader, and then the concentration of BALP can be calculated.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.