

Human beta-ENaC (Amiloride-sensitive sodium channel subunit beta) ELISA Kit | beta-ENaC elisa kit

Human beta-ENaC (Amiloride-sensitive sodium channel subunit beta) ELISA Kit

CATALOG #
AAA249902

REACTIVITY
Human


PRODUCT SYNONYMS

beta-ENaC (Amiloride-sensitive sodium channel subunit beta); Human beta-ENaC (Amiloride-sensitive sodium channel subunit beta) ELISA Kit; Amiloride-sensitive sodium channel subunit beta; Beta-NaCH; Epithelial Na(+) channel subunit beta; Beta-ENaC; ENaCB; Nonvoltage-gated sodium channel 1 subunit beta; SCNEB; SCNN1B; beta-ENaC; beta-ENaC elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of beta-ENaC. No significant cross-reactivity or interference between beta-ENaC and analogues was observed.
Assay Type	Sandwich
Samples	Serum, plasma, tissue homogenates and other biological fluids
Detection Range	0.156-10ng/ml
Sensitivity	0.094ng/ml
Intra-assay Precision	CV<8%
Inter-assay Precision	CV<10%
Preparation and Storage	Store entire kit at 2-8C for short-term. For longer-term, please store the microplate & standard at -20C, while the remaining reagents can be stored at 2-8C

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for beta-ENaC elisa kit	Principle of the Assay: This kit was based on sandwich enzyme-linked immune-sorbent assay technology. Capture antibody was pre-coated onto 96- well plates. And the biotin conjugated antibody was used as detection antibodies. The standards, test samples and biotin conjugated detection antibody were added to the wells subsequently, and washed with wash buffer. HRP-Streptavidin was added and unbound conjugates were washed away with wash buffer. TMB substrates were used to visualize HRP enzymatic reaction. TMB was catalyzed by HRP to produce a blue color product that changed into yellow after adding acidic stop solution. The density of yellow is proportional to the target amount of sample captured in plate. Read the O.D. absorbance at 450nm in a microplate reader, and then the concentration of target can be calculated.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	124301196
NCBI GeneID	6338
NCBI Accession #	NP_000327.2
NCBI GenBank Nucleotide #	NM_000336.2
UniProt Accession #	P51168; O60891; Q96KG2; Q9UJ32; Q9UMU5; C5HTZ2
Molecular Weight	72,659 Da
NCBI Official Full Name	amiloride-sensitive sodium channel subunit beta
NCBI Official Synonym Full Names	sodium channel, non-voltage-gated 1, beta subunit
NCBI Official Symbol	SCNN1B
NCBI Official Synonym Symbols	BESC1; ENaCb; SCNEB; ENaCbeta
NCBI Protein Information	amiloride-sensitive sodium channel subunit beta; beta-ENaC; beta-NaCH; epithelial Na(+) channel subunit beta; epithelial sodium channel beta-2 subunit; epithelial sodium channel beta-3 subunit; sodium channel, nonvoltage-gated 1, beta; nasal epithelial so
UniProt Protein Name	Amiloride-sensitive sodium channel subunit beta
UniProt Gene Name	SCNN1B
UniProt Synonym Gene Names	Beta-ENaC; ENaCB
UniProt Entry Name	SCNNB_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/beta-enac-amiloride-sensitive-sodium-channel-subunit-beta/249902>