

Human IER3 (Immediate Early Response 3) ELISA Kit | IER3 elisa kit

Human IER3 (Immediate Early Response 3) QuickTest ELISA Kit

CATALOG #
AAA252366

REACTIVITY
Human


PRODUCT SYNONYMS

IER3 (Immediate Early Response 3); Human IER3 (Immediate Early Response 3) QuickTest ELISA Kit; Radiation-inducible immediate-early gene IEX-1; Differentiation-dependent gene 2 protein; Protein DIF-2; Immediate early protein GLY96; Immediate early response 3 protein; PACAP-responsive gene 1 protein; Protein PRG1; IER3; DIF2; IEX1; PRG1; IER3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Assay Type	Sandwich ELISA, Double Antibody
Samples	Serum, plasma, cell culture supernatant, cell lysate or tissue lysate, other biological fluid samples
Detection Range	78.125-5000pg/ml
Sensitivity	46.875pg/ml
Time	120 minutes
Preparation and Storage	Store at 2-8 degree C for 12 months

IMAGES



Standard Curve (Sample)

NCBI AND UNIPROT INFORMATION

NCBI GI #	119964723
NCBI GeneID	8870
NCBI Accession #	NP_003888.2
NCBI GenBank Nucleotide #	NM_003897.3
UniProt Accession #	P46695 ; Q5SU30 ; Q92691 ; Q93044
Molecular Weight	16,903 Da
NCBI Official Full Name	radiation-inducible immediate-early gene IEX-1

NCBI Official Synonym Full Names	immediate early response 3
NCBI Official Symbol	IER3
NCBI Official Synonym Symbols	DIF2; IEX1; PRG1; DIF-2; GLY96; IEX-1; IEX-1L
NCBI Protein Information	radiation-inducible immediate-early gene IEX-1; PACAP-responsive gene 1 protein; anti-death protein; differentiation-dependent gene 2 protein; expressed in pancreatic carcinoma; gly96, mouse, homolog of; immediate early protein GLY96; immediate early resp
UniProt Protein Name	Radiation-inducible immediate-early gene IEX-1
UniProt Gene Name	IER3
UniProt Synonym Gene Names	DIF2; IEX1; PRG1; Protein DIF-2; Protein PRG1
UniProt Entry Name	IEX1_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/ier3-immediate-early-response-3/252366>