

General D-Amino Acid Oxidase Assay Kit

D-Amino Acid Oxidase Activity Microplate Assay Kit

CATALOG # AAA266480	REACTIVITY General	APPLICATIONS Functional Assay
-------------------------------	------------------------------	---


PRODUCT SYNONYMS

D-Amino Acid Oxidase; D-Amino Acid Oxidase Activity Microplate Assay Kit; EC 1.4.3.3; DAAO; DAMOX; DAO; D-Amino Acid Oxidase assay kit

PRODUCT SPECIFICATIONS


Reactivity	General
Samples	Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Detection Method	Colorimetric
Detection Range	3umol/L-300umol/L
CAS Number	9000-88-8
Platform	Microplate reader 450nm
Applicable Applications for D-Amino Acid Oxidase assay kit	FA (Functional Assay)
Preparation and Storage	Shipped and store at 4 degree C for 6 months.

IMAGES



Application Data

Positive Control reaction in 96-well plate assay with decreasing the concentration



Standard Curve (Sample)

Representative Standard Curve

ADDITIONAL INFORMATION

Related Product Information for D-Amino Acid Oxidase assay kit	Detection and Quantification of D-Amino Acid Oxidase Activity.
--	--

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.