

Mouse CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit | CDKN2A elisa kit

Mouse CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit

CATALOG #
AAA270138

REACTIVITY
Mouse


PRODUCT SYNONYMS

CDKN2A (Cyclin Dependent Kinase Inhibitor 2A); Mouse CDKN2A (Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit; P16; ARF; CDK4I; CMM2; INK4; INK4a; MLM; MTS1; TP16; P14; P14ARF; P16INK4a; P19; Melanoma P16; Cyclin-dependent kinase 4 inhibitor A; Multiple tumor suppressor 1; CDKN2A elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Samples	Tissue homogenates and other biological fluids.
Assay Type	Sandwich
Detection Range	0.32-20 ng/mL
Sensitivity	0.108 ng/mL
Intra-assay Precision	CV%<8%
Inter-assay Precision	CV%<10%
Preparation and Storage	Store at 2-8 degree C. Stable for 6 months from date of receipt.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CDKN2A elisa kit

Principle of the assay: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Cyclin Dependent Kinase Inhibitor 2A(CDKN2A). Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Cyclin Dependent Kinase Inhibitor 2A(CDKN2A). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Cyclin Dependent Kinase Inhibitor 2A(CDKN2A), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of Cyclin Dependent Kinase Inhibitor 2A(CDKN2A) in the samples is then determined by comparing the OD of the samples to the standard curve.

Product Categories/Family for CDKN2A elisa kit

[Tumor immunity](#); [Developmental science](#)

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/mouse-elisa-kits/cdkn2a-cyclin-dependent-kinase-inhibitor-2a/270138>