

Mouse GIP (Gastric Inhibitory Polypeptide) ELISA Kit | GIP elisa kit

Mouse GIP (Gastric Inhibitory Polypeptide) ELISA Kit

CATALOG #
AAA270237

REACTIVITY
Mouse


PRODUCT SYNONYMS

GIP (Gastric Inhibitory Polypeptide); Mouse GIP (Gastric Inhibitory Polypeptide) ELISA Kit; Incretin hormone; Glucose Dependent Insulinotropic Peptide; GIP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Samples	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Assay Type	Competitive
Detection Range	78.13-5000 pg/mL
Sensitivity	25.2 pg/mL
Intra-assay Precision	CV%<8%
Inter-assay Precision	CV%<10%
Preparation and Storage	Store at 2-8 degree C. Stable for 6 months from date of receipt.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for GIP elisa kit

Principle of the assay: This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with Gastric Inhibitory Polypeptide(GIP) protein. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Gastric Inhibitory Polypeptide(GIP). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of Gastric Inhibitory Polypeptide(GIP) in the samples is then determined by comparing the OD of the samples to the standard curve.

Product Categories/Family for GIP elisa kit

[Endocrinology](#); [Gastroenterology](#); [Hormone metabolism](#)

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

