

Human ProGRP (Pro-Gastrin Releasing Peptide) ELISA Kit | ProGRP elisa kit

Human ProGRP (Pro-Gastrin Releasing Peptide) ELISA Kit

CATALOG #
AAA270301

REACTIVITY
Human


PRODUCT SYNONYMS

ProGRP (Pro-Gastrin Releasing Peptide); Human ProGRP (Pro-Gastrin Releasing Peptide) ELISA Kit; pro-GRP; preproGRP; pre-proGRP; ProGRP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Assay Type	Competitive
Detection Range	15.63-1000 pg/mL
Sensitivity	4.95 pg/mL
Intra-assay Precision	CV%<8%
Inter-assay Precision	CV%<10%
Preparation and Storage	Store at 2-8 degree C. Stable for 6 months from date of receipt.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for ProGRP elisa kit

Principle of the assay: This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with Pro-Gastrin Releasing Peptide(ProGRP) protein. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Pro-Gastrin Releasing Peptide(ProGRP). Next,Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of Pro-Gastrin Releasing Peptide(ProGRP) in the samples is then determined by comparing the OD of the samples to the standard curve.

Product Categories/Family for ProGRP elisa kit

[Tumor immunity](#); [Endocrinology](#); [Reproductive science](#); [Pulmonology](#); [Hepatology](#); [Gastroenterology](#); [Hormone metabolism](#)

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

