

Human CDH17 (Cadherin 17) ELISA Kit | CDH17 elisa kit

Human CDH17 (Cadherin 17) ELISA Kit

CATALOG #
AAA271155

REACTIVITY
Human


PRODUCT SYNONYMS

CDH17 (Cadherin 17); Human CDH17 (Cadherin 17) ELISA Kit; CDH16; HPT-1; HPT1; LI-Cadherin; Cadherin; Liver Intestine; Intestinal peptide-associated transporter HPT-1; Liver-intestine cadherin; CDH17 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Assay Type	Sandwich
Detection Range	31.25-2000 pg/mL
Sensitivity	14.1 pg/mL
Intra-assay Precision	CV%<8%
Inter-assay Precision	CV%<10%
Preparation and Storage	Store at 2-8 degree C. Stable for 6 months from date of receipt.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CDH17 elisa kit

Principle of the assay: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Cadherin 17(CDH17). Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Cadherin 17(CDH17). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Cadherin 17(CDH17), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Cadherin 17(CDH17) in the samples is then determined by comparing the OD of the samples to the standard curve.

Product Categories/Family for CDH17 elisa kit

CD & Adhesion molecule; Tumor immunity; Infection immunity; Hepatology; Gastroenterology

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/cdh17-cadherin-17/271155>