

Horse cAMP (Cyclic Adenosine Monophosphate) ELISA Kit | cAMP elisa kit

Horse cAMP (Cyclic Adenosine Monophosphate) ELISA Kit

CATALOG #
AAA276273

REACTIVITY
Horse


PRODUCT SYNONYMS

cAMP (Cyclic Adenosine Monophosphate); Horse cAMP (Cyclic Adenosine Monophosphate) ELISA Kit; c-AMP; 3'-5'-Cyclic Adenosine Monophosphate; Adenosine Cyclophosphate; cAMP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Horse
Specificity	This assay has high sensitivity and excellent specificity for detection of Horse cAMP. No significant cross-reactivity or interference between Horse cAMP and analogues was observed.
Assay Type	Competitive
Detection Range	312.5-20000 pg/mL
Sensitivity	92.5 pg/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay)? CV%<8%. Three samples of known concentration were tested twenty times on one plate to assess intra-assay precision.
Inter-assay Precision	Inter-assay Precision (Precision between assays)? CV%<10%. Three samples of known concentration were tested in forty separate assays to assess inter-assay precision.
Preparation and Storage	Store at 2-8 degree C. Stable for 6 months from date of receipt.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for cAMP elisa kit

Principle of the Assay: This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with Cyclic Adenosine Monophosphate(cAMP) protein. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Cyclic Adenosine Monophosphate(cAMP). Next,Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Cyclic Adenosine Monophosphate(cAMP) in the samples is then determined by comparing the OD of the samples to the standard curve.

Product Categories/Family for cAMP elisa kit

[Signal transduction](#); [Metabolic pathway](#)

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/horse-elisa-kits/camp-cyclic-adenosine-monophosphate/276273>