

Human CLU (Clusterin) ELISA Kit | CLU elisa kit

Human CLU (Clusterin) ELISA Kit

CATALOG #
AAA278040

REACTIVITY
Human


PRODUCT SYNONYMS

CLU (Clusterin); Human CLU (Clusterin) ELISA Kit; CLU; TRPM2; KUB1; CLI; APOJ; SGP2; AAG4; APO-J; SP-40; Apolipoprotein J; Complement Lysis Inhibitor; 40; Sulfated Glycoprotein 2; Testosterone-Repressed Prostate Message 2; CLU elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Assay Type	Sandwich
Samples	serum, plasma, urine, saliva, tissue homogenates, cerebrospinal fluid, cell culture supernates and other biological fluids
Detection Range	1.57-100ug/mL
Sensitivity	0.57ug/mL

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CLU elisa kit

Test Principle: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Human CLU. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Human CLU. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Human CLU, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm +/- 10nm. The concentration of Human CLU in the samples is then determined by comparing the OD of the samples to the standard curve.

Product Categories/Family for CLU elisa kit

[Apoptosis](#); [Tumor immunity](#); [Infection immunity](#)

NCBI AND UNIPROT INFORMATION

UniProt Accession #	O75153
UniProt Protein Name	Protein KIAA0664

UniProt Gene Name	KIAA0664
UniProt Entry Name	K0664_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/clu-clusterin/278040>