

Rat LpPLA2 (Phospholipase A2, Lipoprotein Associated) ELISA Kit | LpPLA2 elisa kit

Rat LpPLA2 (Phospholipase A2, Lipoprotein Associated) ELISA Kit

CATALOG #
AAA278464

REACTIVITY
Rat


PRODUCT SYNONYMS

LpPLA2 (Phospholipase A2; Lipoprotein Associated); Rat LpPLA2 (Phospholipase A2; Lipoprotein Associated) ELISA Kit; PLA2G7; PAF-AH; PAFAH; Lp-PLA2; LDL-PLA2; Platelet Activating Factor Acetylhydrolase; Plasma; Phospholipase A2; Group VII; LDL-associated phospholipase A2; LpPLA2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rat
Assay Type	Sandwich
Samples	Serum, plasma, tissue homogenates and other biological fluids.
Detection Range	3.13-200ng/mL
Sensitivity	5.1pg/mL

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for LpPLA2 elisa kit

Test Principle: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Rat LpPLA2. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Rat LpPLA2. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Rat LpPLA2, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm +/- 10nm. The concentration of Rat LpPLA2 in the samples is then determined by comparing the OD of the samples to the standard curve.

Product Categories/Family for LpPLA2 elisa kit

Enzyme & Kinase; Metabolic pathway; Infection immunity

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.