

Dog SELE (Selectin, Endothelium) ELISA Kit | SELE elisa kit

Dog SELE (Selectin, Endothelium) ELISA Kit

CATALOG # AAA278539	REACTIVITY Dog
-------------------------------	--------------------------


PRODUCT SYNONYMS

SELE (Selectin; Endothelium); Dog SELE (Selectin; Endothelium) ELISA Kit; CD62E; CD62-E; E-Selectin; ELAM; ELAM1; ESEL; SEL-E; E-LAM; E-LAM1; E-SEL; LECAM2; Endothelial Leukocyte Adhesion Molecule 1; CD62 Antigen-Like Family Member E; SELE elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Dog
Assay Type	Sandwich
Samples	serum, plasma, tissue homogenates and other biological fluids
Detection Range	1.57-100ng/mL
Sensitivity	0.61ng/mL

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for SELE elisa kit	Test Principle: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Dog SELE. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Dog SELE. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Dog SELE, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm +/- 10nm. The concentration of Dog SELE in the samples is then determined by comparing the OD of the samples to the standard curve.
Product Categories/Family for SELE elisa kit	Signal transduction; CD & Adhesion molecule; Tumor immunity; Infection immunity

NCBI AND UNIPROT INFORMATION

UniProt Accession #	P33730
UniProt Protein Name	E-selectin

UniProt Gene Name	SELE
UniProt Synonym Gene Names	ELAM-1; LECAM2

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/dog-elisa-kits/sele-selectin-endothelium/278539>