

Dog SAA (Serum Amyloid A) ELISA Kit | SAA elisa kit

Dog SAA (Serum Amyloid A) ELISA Kit

CATALOG #
AAA278783

REACTIVITY
Dog


PRODUCT SYNONYMS

SAA (Serum Amyloid A); Dog SAA (Serum Amyloid A) ELISA Kit; SA-A; SAA1; PIG4; TP53I4; Amyloid protein A; Amyloid fibril protein AA; Serum amyloid A-1 protein; SAA elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Dog
Assay Type	Sandwich
Samples	serum, plasma and other biological fluids
Detection Range	3.13-200ng/mL
Sensitivity	0.273ng/mL

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for SAA elisa kit

Test Principle: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Dog SAA. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Dog SAA. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Dog SAA, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm +/- 10nm. The concentration of Dog SAA in the samples is then determined by comparing the OD of the samples to the standard curve.

Product Categories/Family for SAA elisa kit

[Infection immunity](#); [Kidney biomarker](#); [Hepatology](#); [Urology](#)

NCBI AND UNIPROT INFORMATION

NCBI GI #	548887
NCBI GeneID	403585
UniProt Accession #	P19708

Molecular Weight	13.8 kDa
NCBI Official Full Name	Serum amyloid A protein
NCBI Official Symbol	SAA1
NCBI Official Synonym Symbols	SAA
NCBI Protein Information	serum amyloid A1
UniProt Protein Name	Serum amyloid A protein
UniProt Gene Name	SAA1
UniProt Synonym Gene Names	SAA
UniProt Entry Name	SAA_CANLF

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/dog-elisa-kits/saa-serum-amyloid-a/278783>