

Rabbit MARCO Polyclonal Antibody | anti-MARCO antibody

MARCO Polyclonal Antibody

CATALOG # AAA280660	HOST Rabbit	APPLICATIONS Western Blot
-------------------------------	-----------------------	-------------------------------------


PRODUCT SYNONYMS

MARCO; Polyclonal Antibody; MARCO Polyclonal Antibody; anti-MARCO antibody

PRODUCT SPECIFICATIONS

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purity/Purification	Affinity purification
Form/Format	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Sequence Length	520
Applicable Applications for anti-MARCO antibody	WB (Western Blot)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 421-520 of human MARCO (NP_006761.1).
Species Reactivity	Human, Mouse, Rat
Preparation and Storage	Store at -20°C. Avoid freeze / thaw cycles.

IMAGES



WB (Western Blot)

Western blot analysis of extracts of various cell lines, using MARCO antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.

ADDITIONAL INFORMATION

Product Categories/Family for anti-MARCO antibody	Polyclonal
---	------------

NCBI AND UNIPROT INFORMATION

NCBI GI #	5803080
NCBI GeneID	8685
NCBI Accession #	NP_006761.1
NCBI GenBank Nucleotide #	NM_006770.3
UniProt Accession #	Q9UEW3; Q9Y5S3; B4DW79
Molecular Weight	Observed MW:60kDa Calculated MW:43kDa/52kDa
NCBI Official Full Name	macrophage receptor MARCO
NCBI Official Synonym Full Names	macrophage receptor with collagenous structure
NCBI Official Symbol	MARCO
NCBI Official Synonym Symbols	SCARA2
NCBI Protein Information	macrophage receptor MARCO
UniProt Protein Name	Macrophage receptor MARCO
UniProt Gene Name	MARCO
UniProt Synonym Gene Names	SCARA2
UniProt Entry Name	MARCO_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.