

Human Luteinizing Hormone (LH) ELISA Kit | LH elisa kit

Human Luteinizing Hormone (LH) ELISA Kit

CATALOG #
AAA282349

REACTIVITY
Human


PRODUCT SYNONYMS

Luteinizing Hormone (LH); Human Luteinizing Hormone (LH) ELISA Kit; LH elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Assay Type	Sandwich ELISA
Sample Type	Serum, plasma and cell culture supernates.
Sensitivity	31pg/mL
Standard Curve Range	62.5-4000pg/mL
Detection	100ul
Assay Time	4hours
Instrument	Microplate Reader
Conjugation	Biotin
Components	Microtiter Plate, Enzyme Conjugate, Standards, Substrate A, Substrate B, Stop Solution, Wash Solution, Balance Solution
Intended Use	This Human Luteinizing Hormone (LH) ELISA kit is intended for laboratory research use only. The stop solution changes color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of Human Luteinizing Hormone (LH) in the sample, this Human Luteinizing Hormone (LH) ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of optical density versus Human Luteinizing Hormone (LH) concentration. The concentration of the samples is then determined by comparing the O.D. of the samples to the standard curve.
Preparation and Storage	Store unopened kit at 2-8 degree C for 1 week. If more than one week, please keep the components of the kit according to the instructions. Do not use past kit expiration date.

IMAGES



Application Data

ADDITIONAL INFORMATION

Related Product Information for LH elisa kit	The ELISA (Enzyme-Linked Immunosorbent Assay) kit is an in vitro enzyme-linked immunosorbent assay for the quantitative measurement of samples in Serum, plasma and cell culture supernates.
--	--

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/luteinizing-hormone-lh/282349>