

Human CD69 ELISA Kit | CD69 elisa kit

Human CD69 ELISA Kit

CATALOG #
AAA282457

REACTIVITY
Human


PRODUCT SYNONYMS

CD69; Human CD69 ELISA Kit; CD69 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	The method has high sensitivity and specificity for the detection of CD69. there was no significant cross-reactivity or interference between CD69 and analogues.
Assay Type	Quantitative Sandwich
Detection Range	46.87-3000pg/mL
Sensitivity	< 23.4pg/mL
Intra-assay Precision	Intra-plate Precision 3 samples with low, middle and high level CD69 were tested 20 times on one plate, respectively. Intra-Assay: CV<10%
Inter-assay Precision	Inter-plate Precision 3 samples with low, middle and high level CD69 were tested on 3 different plates, 20 replicates in each plate. Inter-Assay: CV<15%
Preparation and Storage	Ship: Ice Pack

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CD69 elisa kit	Background: The kit is a sandwich enzyme immunoassay for in vitro quantitative measurement of CD69 in human serum, plasma, cell culture supernatants and other biological fluids. Principle of the Assay: This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for human CD69 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any CD69 present is bound by the immobilized antibody. After washing away any unbound substances, and then a detection antibody specific for CD69 is added to the wells and binds to the combination of capture antibody CD69 in sample. Following a wash to remove any unbound combination, and enzyme conjugate is added to the wells. Following incubation and wash steps, a substrate solution is added to the wells and color develops in proportion to the amount of CD69 bound in the initial step. The color development is stopped and the absorbance is measured.
--	---

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/cd69/282457>