

Human sonic hedgehog N-Terminus, Shh-N ELISA Kit | Shh-N elisa kit

Human sonic hedgehog N-Terminus, Shh-N ELISA Kit

CATALOG #
AAA29406

REACTIVITY
Human

PRODUCT SYNONYMS

sonic hedgehog N-Terminus; Shh-N; Human sonic hedgehog N-Terminus; Shh-N ELISA Kit; Shh-N elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	No significant cross-reactivity or interference between this analyte and analogues is observed.
Assay Type	Quantitative Sandwich
Detection Range	3.12 ng/ml - 100 ng/ml
Sensitivity	1.0 ng/ml.
Intra-assay Precision	Intra-assay CV (%) is less than 15%
Inter-assay Precision	Inter-assay CV (%) is less than 15%. [CV(%) = SD/mean ×100].
Preparation and Storage	Store all reagents at 2-8 degree C

ADDITIONAL INFORMATION

Related Product Information for Shh-N elisa kit

Background: This Quantitative Sandwich ELISA kit is only for in vitro research use only! This kit is intended to be used for determination the level of SHHN (hereafter termed "analyte") in undiluted original Human body fluids, tissue homogenates, secretions or feces samples. This kit is NOT suitable for assaying non-biological sources of substances.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/sonic-hedgehog-n-terminus-shh-n/29406>