

Dog carcinoembryonic antigen, CEA ELISA Kit | CEA elisa kit

Dog Carcinoembryonic antigen, CEA ELISA Kit

CATALOG #
AAA29440

REACTIVITY
Dog

PRODUCT SYNONYMS

carcinoembryonic antigen; CEA; Dog Carcinoembryonic antigen; CEA ELISA Kit; Caninecarcinoembryonic antigen; carcinoembryonic antigen; CEA elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Dog
Assay Type	Quantitative Sandwich
Detection Range	0.5ng/ml-16ng/ml
Sensitivity	0.1ng/ml
Intra-assay Precision	Intra-assay CV (%) is less than 15%.
Inter-assay Precision	Inter-assay CV (%) is less than 15%. [CV(%) = SD/mean ×100]. All CV% should be compared by concentration, not compared by OD values.
Preparation and Storage	Store all reagents at 2-8 degree C

ADDITIONAL INFORMATION

Related Product Information for CEA elisa kit

Background: This Quantitative Sandwich ELISA kit is for lab reagent/research use only, not for drug, household, therapeutic or diagnostic applications! This kit is intended to be used for determine the level of CEAg (hereafter termed "analyte") in undiluted original Dog serum, plasma or tissue homogenates samples. For other sample types please contact tech support to determine compatibility with this assay. This kit is not suitable for assaying non-biological sources of substances.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/dog-elisa-kits/carcinoembryonic-antigen-cea/29440>