

Human KIAA1462 ELISA Kit | KIAA1462 elisa kit

Human Uncharacterized protein KIAA1462, KIAA1462 ELISA Kit

CATALOG #
AAA29480

REACTIVITY
Human

PRODUCT SYNONYMS

KIAA1462; Human Uncharacterized protein KIAA1462; KIAA1462 ELISA Kit; Human Uncharacterized protein KIAA1462 (KIAA1462) ELISA kit; OTTHUMP00000019390; hypothetical protein LOC57608; KIAA1462 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Assay Type	Quantitative Sandwich
Detection Range	0.625ng/ml-20ng/ml
Sensitivity	0.1ng/ml
Intra-assay Precision	Intra-assay CV (%) is less than 15%.
Inter-assay Precision	Inter-assay CV (%) is less than 15%. [CV(%) = SD/mean ×100]. All CV% should be compared by concentration, not compared by OD values.
Preparation and Storage	Store all reagents at 2-8 degree C

ADDITIONAL INFORMATION

Related Product Information for KIAA1462 elisa kit	Background: This Quantitative Sandwich ELISA kit is for lab reagent/research use only, not for drug, household, therapeutic or diagnostic applications! This kit is intended to be used for determine the level of KIAA1462 (hereafter termed "analyte") in undiluted original Human serum, plasma or tissue homogenates samples. For other sample types please contact tech support to determine compatibility with this assay. This kit is not suitable for assaying non-biological sources of substances.
--	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	148806881
NCBI GeneID	57608
NCBI Accession #	NP_065899.1
NCBI GenBank Nucleotide #	NM_020848.2
UniProt Accession #	Q9P266 ; Q5HYA7 ; Q5T992 ; Q86WZ9 ; Q9BYJ2
Molecular Weight	148,351 Da
NCBI Official Full Name	junctional protein associated with coronary artery disease
NCBI Official Synonym Full Names	KIAA1462
NCBI Official Symbol	KIAA1462
NCBI Official Synonym Symbols	JCAD
NCBI Protein Information	junctional protein associated with coronary artery disease

UniProt Protein Name	Junctional protein associated with coronary artery disease
UniProt Gene Name	KIAA1462
UniProt Synonym Gene Names	JCAD; JCAD
UniProt Entry Name	JCAD_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/kiaa1462/29480>