

Human KIAA0317 ELISA Kit | KIAA0317 elisa kit

Human Protein KIAA0317, KIAA0317 ELISA Kit

CATALOG #
AAA29540

REACTIVITY
Human

PRODUCT SYNONYMS

KIAA0317; Human Protein KIAA0317; KIAA0317 ELISA Kit; Human Protein KIAA0317 (KIAA0317) ELISA kit; hypothetical protein LOC9870; KIAA0317 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Assay Type	Quantitative Sandwich
Detection Range	3.12ng/ml-100ng/ml
Sensitivity	1.0ng/ml
Intra-assay Precision	Intra-assay CV (%) is less than 15%.
Inter-assay Precision	Inter-assay CV (%) is less than 15%. [CV(%) = SD/mean ×100]. All CV% should be compared by concentration, not compared by OD values.
Preparation and Storage	Store all reagents at 2-8 degree C

ADDITIONAL INFORMATION

Related Product Information for KIAA0317 elisa kit

Background: This Quantitative Sandwich ELISA kit is for lab reagent/research use only, not for drug, household, therapeutic or diagnostic applications! This kit is intended to be used for determine the level of KIAA0317 (hereafter termed "analyte") in undiluted original Human serum, plasma or tissue homogenates samples. For other sample types please contact tech support to determine compatibility with this assay. This kit is not suitable for assaying non-biological sources of substances.

NCBI AND UNIPROT INFORMATION

NCBI GI #	87116668
NCBI GeneID	9870
NCBI Accession #	NP_001034568.1
NCBI GenBank Nucleotide #	NM_001039479.1
UniProt Accession #	O15033 ; Q7LDY1 ; Q8IYY9 ; B4E2C7
Molecular Weight	94,223 Da
NCBI Official Full Name	apoptosis-resistant E3 ubiquitin protein ligase 1
NCBI Official Synonym Full Names	apoptosis resistant E3 ubiquitin protein ligase 1
NCBI Official Symbol	AREL1
NCBI Official Synonym Symbols	KIAA0317
NCBI Protein Information	apoptosis-resistant E3 ubiquitin protein ligase 1
UniProt Protein Name	Apoptosis-resistant E3 ubiquitin protein ligase 1

UniProt Gene Name	AREL1
UniProt Synonym Gene Names	KIAA0317
UniProt Entry Name	AREL1_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/kiaa0317/29540>