

Human ring finger protein 145 ELISA Kit | RNF145 elisa kit

Human RING finger protein 145, RNF145 ELISA Kit

CATALOG # AAA29555	REACTIVITY Human
------------------------------	----------------------------

PRODUCT SYNONYMS

ring finger protein 145; Human RING finger protein 145; RNF145 ELISA Kit; Human RING finger protein 145 (RNF145) ELISA kit; DKFZp686M11215; FLJ31951; ; RNF145 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Human Serum, Plasma Or Tissue Homogenates Samples
Assay Type	Quantitative Sandwich
Detection Range	0.625ng/ml-20ng/ml
Sensitivity	0.1ng/ml.
Intra-assay Precision	Intra-assay CV (%) is less than 15%
Inter-assay Precision	Inter-assay CV (%) is less than 15%. [CV(%) = SD/mean ×100]. All CV% should be compared by concentration, not compared by OD values.
Preparation and Storage	Store all reagents at 2-8 degree C

ADDITIONAL INFORMATION

Related Product Information for RNF145 elisa kit	Background/Introduction: This Quantitative Sandwich ELISA kit is for lab reagent/research use only! This kit is intended to be used for determine the level of RNF123 (hereafter termed "analyte") in undiluted original Human serum, plasma or tissue homogenates samples. For other sample types please contact tech support to determine compatibility with this assay. This kit is not suitable for assaying non-biological sources of substances.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	313661390
NCBI GeneID	153830
NCBI Accession #	NP_001186309.1
NCBI GenBank Nucleotide #	NM_001199380.1
UniProt Accession #	Q96MT1; Q8IVP7; B7Z903; B7Z949; E7EVI7
Molecular Weight	75,617 Da
NCBI Official Full Name	RING finger protein 145 isoform 1
NCBI Official Synonym Full Names	ring finger protein 145
NCBI Official Symbol	RNF145
NCBI Protein Information	RING finger protein 145
UniProt Protein Name	RING finger protein 145

UniProt Gene Name	RNF145
UniProt Entry Name	RN145_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/rnf145-human-elisa-kits/ring-finger-protein-145/29555>