

Canine Soluble Alpha-Klotho (SAKL) ELISA Kit | SAKL elisa kit

Canine Soluble Alpha-Klotho (SAKL) ELISA Kit

CATALOG #
AAA29607

REACTIVITY
Canine

PRODUCT SYNONYMS

Soluble Alpha-Klotho (SAKL); Canine Soluble Alpha-Klotho (SAKL) ELISA Kit; SAKL elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Assay Type	Quantitative Sandwich
Detection Range	31.2pg/ml-1000pg/ml
Sensitivity	5.0pg/ml
Intra-assay Precision	Intra-assay CV (%) is less than 15%.
Inter-assay Precision	Inter-assay CV (%) is less than 15%. [CV(%) = SD/mean ×100]. All CV% should be compared by concentration, not compared by OD values.
Preparation and Storage	Store all reagents at 2-8 degree C

ADDITIONAL INFORMATION

Related Product Information for SAKL elisa kit

Background: This Quantitative Sandwich ELISA kit is for lab reagent/research use only, not for drug, household, therapeutic or diagnostic applications! This kit is intended to be used for determine the level of SAKL (hereafter termed "analyte") in undiluted original Canine serum, plasma or tissue homogenates samples. For other sample types please contact tech support to determine compatibility with this assay. This kit is not suitable for assaying non-biological sources of substances.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/canine-elisa-kits/soluble-alpha-klotho-sakl/29607>