

Human cyclin G associated kinase ELISA Kit | GAK elisa kit

Human Cyclin-G-associated kinase, GAK ELISA Kit

CATALOG #
AAA298437

REACTIVITY
Human

PRODUCT SYNONYMS

cyclin G associated kinase; Human Cyclin-G-associated kinase; GAK ELISA Kit; Human Cyclin-G-associated kinase (GAK) ELISA kit; FLJ16629; FLJ40395; MGC99654; ; GAK elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	No significant cross-reactivity or interference between this analyte and analogues is observed.
Samples	Body fluids, tissue homogenates, secretions or feces samples
Assay Type	Quantitative Sandwich
Detection Range	0.25 ng/ml - 8 ng/ml
Sensitivity	0.1 ng/ml.
Intra-assay Precision	Intra-assay CV (%) is less than 15%
Inter-assay Precision	Inter-assay CV (%) is less than 15%. [CV(%) = SD/mean ×100].
Preparation and Storage	Store all reagents at 2-8 degree C

ADDITIONAL INFORMATION

Related Product Information for GAK elisa kit	Background/Introduction: This Quantitative Sandwich ELISA kit is only for in vitro research use only! This kit is intended to be used for determination the level of GAK (hereafter termed "analyte") in undiluted original Human body fluids, tissue homogenates, secretions or feces samples. This kit is NOT suitable for assaying non-biological sources of substances.
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	558472758
NCBI GeneID	2580
NCBI Accession #	NP_001273762.1
NCBI GenBank Nucleotide #	NM_001286833.1
UniProt Accession #	O14976; Q5U4P5; Q9BVY6
Molecular Weight	143,191 Da
NCBI Official Full Name	cyclin-G-associated kinase isoform 2
NCBI Official Synonym Full Names	cyclin G associated kinase
NCBI Official Symbol	GAK
NCBI Official Synonym Symbols	DNAJ26; DNAJC26

NCBI Protein Information	cyclin-G-associated kinase; auxilin-2
UniProt Protein Name	Cyclin-G-associated kinase
UniProt Gene Name	GAK
UniProt Entry Name	GAK_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/gak-human-elisa-kits/cyclin-g-associated-kinase/298437>