

Human homeobox A5 ELISA Kit | HOXA5 elisa kit

Human Homeobox protein Hox-A5, HOXA5 ELISA Kit

CATALOG #
AAA298900

REACTIVITY
Human

PRODUCT SYNONYMS

homeobox A5; Human Homeobox protein Hox-A5; HOXA5 ELISA Kit; Human Homeobox protein Hox-A5 (HOXA5) ELISA kit; HOX1; HOX1.3; HOX1C; MGC9376; homeo box 1C; homeo box A5; homeobox protein HOXA5; HOXA5 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	No significant cross-reactivity or interference between this analyte and analogues is observed.
Samples	Rat Body Fluids And Tissue Homogenates
Assay Type	Quantitative Sandwich
Detection Range	0.25ng/ml-8ng/ml
Sensitivity	0.1ng/ml.
Intra-assay Precision	Intra-assay CV (%) is less than 15%
Inter-assay Precision	Inter-assay CV (%) is less than 15%. [CV(%) = SD/mean ×100]. All CV% should be compared by concentration, not compared by OD values.
Preparation and Storage	Store all reagents at 2-8 degree C

ADDITIONAL INFORMATION

Related Product Information for HOXA5 elisa kit	Background/Introduction: This Quantitative Sandwich ELISA kit is for lab reagent/research use only! This kit is intended to be used for determine the level of HO2 (hereafter termed "analyte") in undiluted original Rat body fluids and tissue homogenates. This kit is NOT suitable for assaying non-biological sources of substances.
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	24497517
NCBI GeneID	3202
NCBI Accession #	NP_061975.2
NCBI GenBank Nucleotide #	NM_019102.3
UniProt Accession #	P20719 ; O43367 ; Q96CY6 ; A4D179
Molecular Weight	29,345 Da
NCBI Official Full Name	homeobox protein Hox-A5
NCBI Official Synonym Full Names	homeobox A5
NCBI Official Symbol	HOXA5
NCBI Official Synonym Symbols	HOX1; HOX1C; HOX1.3

NCBI Protein Information	homeobox protein Hox-A5; homeo box 1C; homeo box A5; homeobox protein HOXA5; homeobox protein Hox-1C
UniProt Protein Name	Homeobox protein Hox-A5
UniProt Gene Name	HOXA5
UniProt Synonym Gene Names	HOX1C
UniProt Entry Name	HXA5_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/hoxa5-human-elisa-kits/homeobox-a5/298900>