

Human family with sequence similarity 161, member A ELISA Kit | FAM161A elisa kit

Human Protein FAM161A, FAM161A ELISA Kit

CATALOG #
AAA299307

REACTIVITY
Human

PRODUCT SYNONYMS

family with sequence similarity 161; member A; Human Protein FAM161A; FAM161A ELISA Kit; Human Protein FAM161A (FAM161A) ELISA kit; FLJ13305; MGC129982; MGC129983; hypothetical protein LOC84140; FAM161A elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Undiluted original Human serum, plasma or tissue homogenates samples
Assay Type	Quantitative Sandwich
Detection Range	0.25ng/ml-8ng/ml
Sensitivity	0.1ng/ml
Intra-assay Precision	Intra-assay CV (%) is less than 15%.
Inter-assay Precision	Inter-assay CV (%) is less than 15%. [CV(%) = SD/mean ×100]. All CV% should be compared by concentration, not compared by OD values.
Preparation and Storage	Store all reagents at 2-8 degree C

ADDITIONAL INFORMATION

Related Product Information for FAM161A elisa kit	Background/Introduction: This Quantitative Sandwich ELISA kit is for lab reagent/research use only! This kit is intended to be used for determine the level of FAM161A (hereafter termed "analyte") in undiluted original Human serum, plasma or tissue homogenates samples. For other sample types please contact tech support to determine compatibility with this assay. This kit is not suitable for assaying non-biological sources of substances.
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	320089540
NCBI GeneID	84140
NCBI Accession #	NP_001188472.1
NCBI GenBank Nucleotide #	NM_001201543.1
UniProt Accession #	Q3B820; Q9H8R2; B4DJV7
Molecular Weight	76,752 Da
NCBI Official Full Name	protein FAM161A isoform 1
NCBI Official Synonym Full Names	family with sequence similarity 161, member A
NCBI Official Symbol	FAM161A
NCBI Official Synonym Symbols	RP28
NCBI Protein Information	protein FAM161A; retinitis pigmentosa 28 (autosomal recessive)

UniProt Protein Name	Protein FAM161A
UniProt Gene Name	FAM161A
UniProt Entry Name	F161A_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/fam161a-human-elisa-kits/family-with-sequence-similarity-161-member-a/299307>