

Nucleoprotein Recombinant Protein | N recombinant protein

Recombinant Crimean-Congo hemorrhagic fever virus Nucleoprotein

CATALOG # AAA309804	HOST E Coli	REACTIVITY Crimean-Congo hemorrhagic fever virus
-------------------------------	-----------------------	--

PRODUCT SYNONYMS

Nucleoprotein; Recombinant Crimean-Congo hemorrhagic fever virus Nucleoprotein; Nucleocapsid protein; N recombinant protein

PRODUCT SPECIFICATIONS

Host	E Coli
Reactivity	Crimean-Congo hemorrhagic fever virus
Purity/Purification	Greater than 90% as determined by SDS-PAGE.
Form/Format	Tris-based buffer 50% glycerol
Sequence Positions	Full Length, 1-482aa
Sequence	MENKIEVNNK DEMNKFEEF KKGNGLVDTF TNPYSFCESV PNLERFVFQM ASATDDAQKD SIYASALVEA TKFCAPIYEC AWWSSSTGIVK KGLEWFEKNA GTIKSWDESY IELKVEVPKI EQLANYQQAA LKWRKDIGFR VNANTAALSH KVLAEYKVPK EIVMSVKEML SDMIRRRNLI LNRGGDENPR GPVSREHVEW CREFVKGKYI MAFNPPWGD I NKSGRSGIAL VATGLAKLAE TEGKGVFDEA KKTVEALNGY LDKHKDEV DK ASADNMITNL LKHIKAQEL YKNSSALRAQ GAQIDTAFSS YYWLYKAGVT PETFPTVSQF LFELGKQPRG TKKMKKALLS TPMKWGKKLY ELFADDSFQQ NRIYMHPAVL TAGRISEMGV CFGTIPVANP DDAAQGS GHT KSILNLRNT ETNNPCA KTI VKLFEIQKTG FNIQDMDIVA SEHLLHQSLV GKQSPFQ NAY NVKGNATSAN II
Immunogen Description	Expression Region:1-482 aa; Sequence Info:Full Length
Calculated MW	68 kDa.
Tag Info	N-terminal 6xHis-B2M-tagged
Preparation and Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 monthsat -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week

ADDITIONAL INFORMATION

Related Product Information for N recombinant protein	RNA binding
---	-------------

NCBI AND UNIPROT INFORMATION

NCBI GI #	127869
UniProt Accession #	P27317
NCBI Official Full Name	Nucleoprotein
UniProt Protein Name	Nucleoprotein
UniProt Gene Name	N

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/n-recombinant-protein-e-coli/309804>