

Human Histone H3 ELISA Kit | H3 elisa kit

Histone H3 (Di Methyl-Lys18) Colorimetric Cell-Based ELISA Kit

CATALOG #
AAA319327

REACTIVITY
Human

PRODUCT SYNONYMS

Histone H3; Histone H3 (Di Methyl-Lys18) Colorimetric Cell-Based ELISA Kit; Histone H3 (Di Methyl-Lys18) Cell-Based ELISA Kit; H3K18ME2; HIST1H3A; H3FA; HIST1H3B; H3FL; HIST1H3C; H3FC; HIST1H3D; H3FB; HIST1H3E; H3FD; HIST1H3F; H3FI; HIST1H3G; H3FH; HIST1H3H; H3FK; HIST1H3I; H3FF; HIST1H3J; H3FJ; Histone H3.1; Histone H3/a; Histone H3/b; Histone H3/c; Histone H3/d; Histone H3/f; Histone H3/h; Histone H3/i; Histone H3/j; Histone H3/k; Histone H3/l; HIST2H3A; HIST2H3C; H3F2; H3FM; HIST2H3D; Histone H3.2; Histone H3/m; Histone H3/o; H3F3A; H3.3A; H3F3; PP781; H3F3B; H3.3B; Histone H3.3; YM3333; HIST1H3A/HIST1H3B/HIST1H3C/HIST1H3D/HIST1H3E/HIST1H3F/HIST1H3G/HIST1H3H/HIST1H3I/HIST1H3J/HIST2H3A/HIST2H3C/HIST2H3D/H3F3A/H3F3B; H3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Form/Format	Two 96-Well Plates
Assay Type	Cell-Based ELISA
Samples	Researchers plate their cell line of choice
Detection Range	At least 5000 cells/well
Detection Method	Colorimetric 450 nm
Preparation and Storage	4 degree C/6 Months

IMAGES



WB (Western Blot)

Western blot analysis of HeLa, diluted at 1:1000. Secondary antibody was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent Biotech, MN, USA).

ADDITIONAL INFORMATION

Related Product Information for H3 elisa kit

The Histone H3 (Di Methyl-Lys18) Antibody Colorimetric Cell-Based ELISA Kit is a convenient, lysate-free, high throughput and sensitive assay kit that can detect Histone H3 (Di Methyl-Lys18) Antibody protein expression profile in cells. The kit can be used for measuring the relative amounts of Histone H3 (Di Methyl-Lys18) Antibody in cultured cells as well as screening for the effects that various treatments, inhibitors (ie siRNA or chemicals), or activators have on Histone H3 (Di Methyl-Lys18) Antibody.

NCBI AND UNIPROT INFORMATION

NCBI GeneID	8350/8351/8352/8353/8354/8355/8356/8357/8358/8968
-------------	---

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/histone-h3/319327>