

Histone H4 ELISA Kit | H4 elisa kit

Histone H4 (Acetyl-Lys16) Colorimetric Cell-Based ELISA Kit

CATALOG #
AAA319413

REACTIVITY
Human, Mouse, Rat


PRODUCT SYNONYMS

Histone H4; Histone H4 (Acetyl-Lys16) Colorimetric Cell-Based ELISA Kit; Histone H4 (Acetyl-Lys16) Cell-Based ELISA Kit; HIST1H4A; H4/A; H4FA; HIST1H4B; H4/I; H4FI; HIST1H4C; H4/G; H4FG; HIST1H4D; H4/B; H4FB; HIST1H4E; H4/J; H4FJ; HIST1H4F; H4/C; H4FC; HIST1H4H; H4/H; H4FH; HIST1H4I; H4/M; H4FM; HIST1H4J; H4/E; H4FE; HIST1H4K; H4/D; H4FD; HIST1H4L; H4/K; H4FK; H4k16AC; HIST1H4A/HIST1H4B/HIST1H4C/HIST1H4D/HIST1H4E/HIST1H4F/HIST1H4H/HIST1H4I/HIST1H4J/HIST1H4K/HIST1H4L/HIST2H4A/HIST2H4B/HIST4H4; H4 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human, Mouse, Rat
Form/Format	Two 96-Well Plates
Assay Type	Cell-Based ELISA
Samples	Researchers plate their cell line of choice
Detection Range	At least 5000 cells/well
Detection Method	Colorimetric 450 nm
Preparation and Storage	4 degree C/6 Months

IMAGES



WB (Western Blot)

Western Blot analysis of various cells using Acetyl-Histone H4 (K16) Polyclonal Antibody.
Secondary antibody was diluted at 1:20000

ADDITIONAL INFORMATION

Related Product Information for H4 elisa kit

The Histone H4 (Acetyl-Lys16) Antibody Colorimetric Cell-Based ELISA Kit is a convenient, lysate-free, high throughput and sensitive assay kit that can detect Histone H4 (Acetyl-Lys16) Antibody protein expression profile in cells. The kit can be used for measuring the relative amounts of Histone H4 (Acetyl-Lys16) Antibody in cultured cells as well as screening for the effects that various treatments, inhibitors (ie siRNA or chemicals), or activators have on Histone H4 (Acetyl-Lys16) Antibody.

NCBI AND UNIPROT INFORMATION

NCBI GeneID

121504/554313/8294/8359/8360/8361/8362/8363/8364/8365/8366/8367/8368/8370

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

