

Tissue Fructose (FT) Assay Kit | FT assay kit

Micro Plant Tissue Fructose (FT) Content Assay Kit

CATALOG #
AAA332635

PRODUCT SYNONYMS

Tissue Fructose (FT); Micro Plant Tissue Fructose (FT) Content Assay Kit; FT assay kit

PRODUCT SPECIFICATIONS

Kit Components	Extraction Buffer
	Reagent I
	Reagent II
	Reagent III
	Reagent IV
Features & Benefits	Provides a simple, sensitive, rapid fructose (FT) detection method in plant tissues.
	Detailed sample preparation and result calculation methods are provided.
Usage Notes	Do not mix the components between different batch numbers and manufacturers; otherwise, the results may be abnormal.
	Avoid bubbles while mixing or redissolving components.
	Change pipette tips frequently to avoid cross contamination between components.
	Ensure that all components and equipment are at the proper temperature before starting the experiment.
Preparation and Storage	Please store the whole kit at 4 degree C upon receipt, and protect from light. The kit is valid for 12 months. Please refer to the table below Materials supplied and Storage conditions to store all the components. It is recommended to use fresh samples, which can be stored at -80 degree C for several weeks if the experiment is not carried out immediately. When measuring, the temperature and time of thawing should be controlled. When thawing at room temperature, the sample should be thawed within 4 h.

ADDITIONAL INFORMATION

Related Product Information for FT assay kit	Plant Tissue Fructose (FT)
Product Categories/Family for FT assay kit	Cell analysis; Cell metabolism detection kit; Cytology; Cell Research Tools; Cell Metabolism

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.