

## Acid Protease (ACP) Assay Kit | ACP assay kit

Micro Acid Protease (ACP) Activity Assay Kit

CATALOG #  
**AAA332690**


### PRODUCT SYNONYMS

Acid Protease (ACP); Micro Acid Protease (ACP) Activity Assay Kit; ACP assay kit

### PRODUCT SPECIFICATIONS

|                                |                                                                                                                                                                                                                                                                                                                                                                                                         |
|--------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Kit Components</b>          | Liquid ?<br>Liquid II<br>Reagent II<br>Reagent III<br>Reagent IV<br>Reagent V<br>Standard                                                                                                                                                                                                                                                                                                               |
| <b>Features &amp; Benefits</b> | The verified samples are complete in type and easy to operate                                                                                                                                                                                                                                                                                                                                           |
| <b>Usage Notes</b>             | Fresh samples are recommended and can be stored at - 80 degree C for several weeks if experiments are not performed immediately. When measuring, the temperature and time of thawing should be controlled. When thawing at room temperature, the sample should be thawed within 4 hours.<br>Ensure that all components and equipment are at the appropriate temperature before starting the experiment. |
| <b>Preparation and Storage</b> | Stored at 4 degree C for 6 months, protected from light<br>Shipping: Gel pack with blue ice.                                                                                                                                                                                                                                                                                                            |

### IMAGES



Application Data

### ADDITIONAL INFORMATION

|                                                      |                                                                                                                                                                                                                                               |
|------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Related Product Information for ACP assay kit</b> | Acidic proteinase is an enzyme that catalyzes the hydrolysis of proteins in acidic environments. The enzyme is mainly used in alcohol fermentation, beer brewing, fur softening, fruit wine clarification, soy sauce brewing, feed and so on. |
| <b>Product Categories/Family for ACP assay kit</b>   | Cell analysis; Cell metabolism detection kit                                                                                                                                                                                                  |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.