

Human Growth Differentiation Factor 5 (GDF5) ELISA Kit | GDF5 elisa kit

Human Growth Differentiation Factor 5 (GDF5) ELISA Kit

CATALOG #
AAA43133

REACTIVITY
Human


PRODUCT SYNONYMS

Growth Differentiation Factor 5 (GDF5); Human Growth Differentiation Factor 5 (GDF5) ELISA Kit; GDF5 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	351
Samples	Serum, plasma, cell culture supernates, cell lysates, tissue homogenates
Assay Type	Quantitative Sandwich
Detection Range	0.05ng/ml - 10ng/ml
Sensitivity	0.026ng/ml
Intra-Assay	CV<8%
Inter-Assay	CV<10%
Preparation and Storage	Store at 2 to 8 degree C for 6 months.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for GDF5 elisa kit

Intended Uses: This sandwich kit is for the accurate quantitative detection of Human Growth Differentiation Factor 5 (also known as GDF5) in serum, plasma, cell culture supernates, cell lysates, tissue homogenates.

Principles of the Assay: This kit is an Enzyme-Linked Immunosorbent Assay (ELISA). The plate has been pre-coated with Human GDF5 antibody. GDF5 present in the sample is added and binds to antibodies coated on the wells. And then biotinylated Human GDF5 Antibody is added and binds to GDF5 in the sample. Then Streptavidin-HRP is added and binds to the Biotinylated GDF5 antibody. After incubation unbound Streptavidin-HRP is washed away during a washing step. Substrate solution is then added and color develops in proportion to the amount of Human GDF5. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

NCBI AND UNIPROT INFORMATION

NCBI GI #	215273634
Molecular Weight	39,932 Da
NCBI Official Full Name	growth differentiation factor 5, partial
UniProt Protein Name	Growth differentiation factor 5
UniProt Gene Name	Gdf5

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/growth-differentiation-factor-5-gdf5/43133>